



Introducing YL2800T

All-In-Ones GC Autosampler



YOUNG IN Chromass

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (812)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Уде (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://youngin.nt-rt.ru/> || y nb@nt-rt.ru

M@2800T (1/2)

Overview



- **High sample capacity**
 - HS: 42 samples
 - Liquid: 121 samples
 - SPME: 42 samples
- YL2800T **modality switch operations** were already easier than its competitors: YL2800T is even better though the new **touch screen wizard**
- Features all goodies already available with YL2000H/YL3000A platform:
 - **Quick Installation**
 - **Target Illumination**
 - **Preventive maintenance counters**
 - **CFR 21 Part 11**
- Share the same mounting kit of YL2000H/YL3000A: **5 minutes is the time required to move M@2800T from one GC to another one's**

M@2800T (2/2)

Overview



- **SPME mode:**
 - Re-engineering of **SPME fiber cleaning module** in order to have it better integrated in the design.
 - Added support to 2cm SPME fiber for usage with low concentration compounds
 - Unified inlet for SPME fiber cleaning and HS syringe purge
 - New DEMO mode with temporary license!!
- **Liquid mode:**
 - From 5 to 6 solvents
 - Include double wash mode
- **Headspace: mode:** includes all improvements announced for YL2000H, such as the lowest total cost of ownership in the market

M@2800T: Smart Technologies For The Best User Experience!



- Our R&D team has put together a **unique portfolio of patented, proprietary or licensed technologies.**
- To provide additional robustness to your **headspace analysis:**
 - system integrity test
 - vial leakage.
 - headspace syringes that provide an excellent performance over a large temperature range
- To provide additional robustness to your **liquid sample analysis:**
 - SyringeID
- To provide additional robustness to your **SPME analysis:**
 - Oven design offers superior thermal stability also on long extraction time method

GC Autosamplers

Headspace and All-In-Ones Autosampler

	M@2100H	M@2000H	M@2800T
Headspace	YES (14 samples)	YES (42 samples)	YES (42 samples)
Liquid injection	-	-	YES (121 samples)
SPME	-	-	Optional (42 samples)
User Interface	Keypad	Touchscreen	Touchscreen
Oven position(s)	1	6	6
Shaking capability	YES (Sussultatory)	YES (Orbital)	YES (Orbital)
Programmable injection volume	YES	YES	YES
Supported Headspace Syringe	1, 2.5, 5mL	1, 2.5, 5mL	1, 2.5 mL

Product Models



M@800T

- The standard model include Headspace and Liquid kit. SPME can be sourced separately, also at a later time.

M@850T

- Same as YL2800T. It also includes SPME cleaning device

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://youngin.nt-rt.ru/> || y nb@nt-rt.ru